

Transport of ventilated patients within the hospital

www.medicars.com P.O. Box 49 3600 AA Maarssen Netherlands Tel. +31 (0)30 2613500



Cart for transport of ventilated patients within the hospital

While moving ventilated patients -for instance between Intensive Care Unit (ICU) and operating room (OR) or X-ray department-, the respiratory equipment and any additional equipment as infusion pumps needs to go along with the patient during transport. This cart is at the bedside of the patient and includes all associated equipment necessary during transport. See the next page for a explanation of the configuration.

Configuration

- ✓ SVR 54-150 base cart with side covers, serum holder and bottle brackets
- ✓ Medical isolation transformer with Jobarco Careflex cable and 7-way socket strip
- ✓ Custom battery pack in undercarriage
- ✓ Customized bracket for connection onto the bed
- ✓ Various custom brackets, runners and holders for placing the equipment

Interested?

Feel free to ask for more information about our products. We will be happy to assist you.

For over 40 years Jansen Medicars is manufacturer of medical furniture, monitor mountings and carts



www.medicars.com P.O. Box 49 3600 AA Maarssen Netherlands Tel. +31 (0)30 2613500

- ✓ Serum holder 2x —
- ✓ Philips IntelliVue MP40/MP50 patient monitor
- ✓ Standard runners (6x) for Syringe -Pumps (i.e. (Asena TM TIVA)
- ✓ Dräger Oxylog 3000 plus ventilator
- ✓ Bed connection system ———
- ✓ Battery pack for energy supply during transport and medical isolation transformer

- ✓ Holder for 5 litre gas bottle _____
- ✓ Weinmann Suction Pump WM10700 -





