

## Compact transport cart

While moving ventilated patients -for instance between Intensive Care Unit (ICU) and operating room (OR) or X-ray department-, the respiratory equipment and any additional equipment as infusion pumps needs to go along with the patient during transport. This cart is at the bedside of the patient and includes all associated equipment necessary during transport. See the next page for a explanation of the configuration.



### Configuration

- Flexx two 54-150 basic ICU cart
- Infusion pole and holder for 3 gas cylinders
- Medical grade power cord (Jobarco Careflex) en 6-way power strip
- Bracket for connection onto the bed
- Various custom brackets, runners and holders for placing the equipment

### Interested?

Please feel free to contact us for information or a quotation. We'll be happy to assist you in finding the most optimum configuration.

+31 (0)30 2613500  
[www.medicars.com](http://www.medicars.com)




*Flexx two configuratie*



**Jansen**  
medicars

# TRANSPORT OF VENTILATED PATIENTS WITHIN THE HOSPITAL

- ✓ Serumholder 2x on swivel arm
- ✓ Patient monitor (on pivot)
- ✓ 6-way power strip with 5m careflex power cord
- ✓ Standard runner (3x) for mounting medical devices (syringe pumps, etc.)
- ✓ Dräger Oxilog 3000 ventilator
- ✓ Bed connection system (Hill-rom totalcare bed & Linet bed)
- ✓ 4x 125mm swivel castors, without breaks



- ✓ 2x cable & hose holder
- ✓ Case for emergency facilities intubation equipment and medication
- ✓ Push bar
- ✓ Holder for 2x 5 litre gas cylinder Ø150mm & 1x 2litre cylinder

